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File Edit View Favorites Tools Help

Back Forward Stop Reload Home

Search Favorites Media

Print Mail

W

Address http://www.cdc.gov/niosh/npptl/topics/respirators/disp\_part/ Go Links

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CDC Home CDC Search CDC Health Topics A-Z

NIOSH National Institute for Occupational Safety and Health

Search NIOSH | NIOSH Home | NIOSH Topics | Site Index | Databases and Information Resources | NIOSH Products | Contact Us

NPPTL The National Personal Protective Technology Laboratory

# NIOSH-Approved Disposable Particulate Respirators (Filtering Facepieces)

Recent CDC infection control guidance documents provide recommendations that healthcare workers protect themselves from diseases potentially spread through the air (such as SARS or tuberculosis) by wearing a fit tested respirator at least as protective as a NIOSH-approved N95 respirator.

An N95 respirator is one of nine types of disposable particulate respirators. Particulate respirators are also known as "air-purifying respirators" because they protect by filtering particles out of the air you breathe. Workers can wear any one of the particulate respirators for protection against diseases spread through the air - if they are NIOSH approved and if they have been properly fit tested and maintained. NIOSH-approved disposable respirators are marked with the manufacturer's name, the part number (P/N), the protection provided by the filter (e.g. N95), and "NIOSH."

Select a type of disposable particulate respirator (filtering facepiece) to view a list of manufacturers/suppliers and model numbers:

Type	Description
<a href="#">N95</a>	Filters at least 95% of airborne particles. Not resistant to oil.
<a href="#">N99</a>	Filters at least 99% of airborne particles. Not resistant to oil.
<a href="#">N100</a>	Filters at least 99.97% of airborne particles. Not resistant to oil.
<a href="#">R95</a>	Filters at least 95% of airborne particles. Somewhat resistant to oil.
<a href="#">R99*</a>	Filters at least 99% of airborne particles. Somewhat resistant to oil.
<a href="#">R100*</a>	Filters at least 99.97% of airborne particles. Somewhat resistant to oil.

NPPTL Directory

- [NPPTL Home](#)
- [About NPPTL](#)
- [Contact NPPTL](#)
- [Employment](#)
- [Announcements](#)
- [NPPTL Site Index](#)

Resources

- CBRN Approved Respirators ([APR/APER](#), [SCBA](#))
- [CBRN Respirator Standards Development](#)
- [Other Respirator Standards Development](#)
- [Certified Equipment List](#)
- [Safety and Health Topic: Respirators](#)
- [Safety and Health Topic:](#)

Done

Internet

start Microsoft PowerPoint ... NIOSH-Approved Dis...

08:13 a.m.

[http://www.cdc.gov/niosh/npptl/topics/respirators/disp\\_part/](http://www.cdc.gov/niosh/npptl/topics/respirators/disp_part/)

## NPPTL The National Personal Protective Technology Laboratory

# NIOSH-Approved Disposable Particulate Respirators (Filtering Facepieces)

## NIOSH-Approved N95 Disposable Particulate Respirators

Check the manufacturer for availability of any listed products. Links to manufacturers' web sites are only provided as a convenience to the user. When a mask is both cleared by the FDA as a surgical mask and certified by NIOSH as a N-95 respirator mask the Model Number/Product Line and Approval Numbers are highlighted below in **RED**. On the FDA website search page the mask cleared by the FDA as a surgical mask and certified by NIOSH as an N-95 respirator mask are referred to as a "Surgical N-95 Respirator"

Supplier/ Manufacturer	Contact Number	Model Number/ Product Line	Approval Number
<a href="#">3M Company</a>	1-888-3M HELPS	<b>1860</b>	<b>84A-0006</b>
<a href="#">3M Company</a>	1-888-3M HELPS	8110S, 8210, 8216, 8217, 7048	84A-0007
<a href="#">3M Company</a>	1-888-3M HELPS	N95	84A-0008
<a href="#">3M Company</a>	1-888-3M HELPS	8212 or 8214 respirator	84A-0454
<a href="#">3M Company</a>	1-888-3M HELPS	8511, 8516, 8517, 8616, 8617, 8211	84A-1299
<a href="#">3M Company</a>	1-888-3M HELPS	9211 N95	84A-2668
<a href="#">3M Company</a>	1-888-3M HELPS	9210 N95	84A-2669
<a href="#">3M Company</a>	1-888-3M HELPS	<b>1870 filtering facepiece</b>	<b>84A-3844</b>

## NPPTL Directory

- ▶ [NPPTL Home](#)
- ▶ [About NPPTL](#)
- ▶ [Contact NPPTL](#)
- ▶ [Employment](#)
- ▶ [Announcements](#)
- ▶ [NPPTL Site Index](#)

## Resources

- ▶ CBRN Approved Respirators ([APR/APER](#), [SCBA](#))
- ▶ [CBRN Respirator Standards Development](#)
- ▶ [Other Respirator Standards Development](#)
- ▶ [Certified Equipment List](#)
- ▶ [Safety and Health Topic: Respirators](#)
- ▶ [Safety and Health Topic:](#)

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