

---

---

---


---

---

---

---

---



---

---

---


---

---

---

---

---



---

---

---

---

---

---

---

---

### Specimen Types

Almost any source is acceptable for AFB culture and smear

- **Sputum**
  - **Induced**
  - **Expectorated**
- Bronchial washing/BAL
- Gastric aspirate
- Fresh tissue
- Bone
- Blood
- Bone Marrow
- CSF
- Body fluids
- Abscess
- Stool
- Urine
- Skin



---

---

---

---

---

---

---

---

### Sputum

- Recently discharged material from the bronchial tree, with minimal amounts of oral or nasal material
- Expectorated: from deep productive cough
- Induced: use of nebulization to increase fluid in the airway and ease clearance of sputum

---

---

---

---

---

---

---

---

### Indications for sputum collection:

- To establish an initial diagnosis of TB
- To monitor the infectiousness of the patient
- To determine the effectiveness of treatment



Image Credit: WHO

---

---

---

---

---

---

---

---

### Specimen collection

- Supervise patient for at least the first specimen, until ability to properly collect the specimen has been demonstrated
- Patient should be in a negative pressure room
- Anyone in the room should wear a fit tested N-95 respirator
- All specimens are collected into sealed leak proof containers
- Label specimen with two patient identifiers, collection date/time and specimen type.

---

---

---

---

---

---

---

---

### Specimen collection

- Optimal: collect a diagnostic specimen before the initiation of drug therapy
  - Collect a series of three sputum specimens, 8-24 hours apart, at least one of which is an early morning specimen
- Monitoring of therapy: Obtain sputum specimens for culture at least monthly until cultures convert to negative for TB

---

---

---

---

---


---

---

---

### Specimen collection kits

WSLH Kit #8: Sputum collection  
Order: 1-800-862-1088  
**KITS ARE FREE!**



- Sterile tube with label
- Absorbent pad
- Specimen transport bag
- Cold pack
- Instruction sheet
- Insulated mailer with labels

---

---

---

---

---

---

---

---

### Storage and Transport

- Sputum samples should be refrigerated if they cannot be transported immediately
- Deliver specimens to the laboratory as soon as possible—try not to batch!
- Recommended: Include a cold pack during specimen transport.

---

---

---

---

---

---

---

---

### Storage and Transport

WHY is this important?

- Minimize overgrowth of normal flora
- Viability of AFB
- Rapid turn-around times
  - Isolation precautions
  - Start/Stop treatment

---

---

---

---

---

---

---

---

### Submission of Specimens to WSLH



- Requisition form A
- Order: 1-800-862-1088
- Preprinted with account information
- One form per specimen

---

---

---

---

---

---

---

---

### Submission of Specimens to WSLH

Code	Test Description
MM00250	Mycobacteria (AFB) Smear and Culture
MM00253	Mycobacteria Isolate Identification
MM02881	Mycobacterium tuberculosis Isolate Genotyping
MM00204	Mycobacterium tuberculosis Susceptibility-1st Line Drugs
MM00202	Mycobacterium avium Complex (MAC) Susceptibility
MM00207	Mycobacteria Rapid Grower Susceptibility Organism ID
MM00260	Mycobacterium avium Complex PCR Decontaminated? Yes No Smear Result
MM00256	Mycobacterium tuberculosis PCR Decontaminated? Yes No Smear Result

---

---

---

---

---

---

---

---

### Courier information

Service is offered at no charge to submitters

Call to set up an account and schedule a pickup:  
763-233-0099



---

---

---

---

---

---

---

---

### Assessing Sputum Quality

**Test results are used as an aid in patient diagnosis and treatment.**

**Test results are directly related to the quality of the specimen.**

---

---

---

---

---

---

---

---

### Assessing Specimen Quality

Collection Date/Time:

- CDC MMWR 2005: 54, RR-17: "Persons requiring sputum collection for smear and culture should have at least three consecutive sputum specimens obtained, each collected in 8-24 hour intervals, with at least one being an early morning specimen"
- Minimize transport time to laboratory (no batching!)
- Specimens >7 days old are rejected

---

---

---

---

---

---

---

---

### Assessing Specimen Quality

Sputum Quality

- Specimens are thick and contain mucopurulent material
- 3-5 ml in volume (ideal), but smaller quantities accepted if the quality is satisfactory
- Poor quality specimens are thin and watery—Saliva and nasal secretions are unacceptable
- Induced sputum should be indicated on requisition form to avoid rejection



---

---

---

---





---

---

---

---

### Assessing Specimen Quality

Thick Mucopurulent		Watery (induced?)	
Hemoptysis		Salivary	

<http://www.fda.gov>

---

---

---

---

---

---

---

---

**Mycobacteriology at WSLH**

- Full-service mycobacteriology laboratory
- BSL-3 facility
- Roles:
  - Primary Diagnostic Facility
  - Reference laboratory
  - Public Health Laboratory



---

---

---

---

---

---

---

---

**Mycobacteriology at WSLH**

- 22 laboratories around the state perform smear and culture
  - 4 labs (other than WSLH) perform some level of identification
  - Most others send to WSLH for identification
- WSLH receives clinical specimens from:
  - 2 large local hospitals
  - Health Departments (state-wide)
    - Milwaukee City TB Clinic
    - Madison Dane County Public Health

---

---

---

---

---

---

---

---

**Mycobacteriology at WSLH**

- Smear Microscopy
- PCR for direct detection (NAAT)
- Culture
- Identification
- Drug Susceptibility Testing
  - Molecular
  - Conventional

---

---

---


---

---

---

---

---



---

---

---

---

---


---

---

---

### Smear Microscopy

- Small amount of concentrated patient specimen is placed on slide and stained with Auramine O fluorescent stain
- Rapid and inexpensive screening tool
  - First indication of mycobacterial infection and possible TB disease
  - Must be accompanied by additional testing including culture for confirmatory diagnosis



A photograph of a glass slide with two circular smears of concentrated patient specimen, used for Auramine O staining.

---

---

---

---


---

---

---

---

### Auramine O Smear



A fluorescence micrograph showing bright, rod-shaped structures against a dark background, characteristic of mycobacteria stained with Auramine O.

Photo Credit: laboratoryinfo.com

---

---

---

---

---

---

---

---

**Smear Microscopy: Result Interpretation**

WSLH Report	Graded Scale	Qualitative Scale	Interpretation
Negative	Negative	Negative	Potentially infectious
Rare (1-9 AFB per 100 fields)	1+	Positive	Low-level infectious
Few (1-9 AFB per 10 fields)	2+	Positive	Moderately infectious
Moderate (1-9 AFB per field)	3+	Positive	
Many (>9 AFB per field)	4+	Positive	Highly infectious

---

---

---

---

---

---

---

---

- Smear Microscopy: Limitations**
- Does not distinguish between viable and dead organisms
  - Limited sensitivity
    - High bacterial load: 5,000-10,000 AFB/mL is required for detection
    - Misses >45% of U.S. TB cases
  - Limited specificity
    - All mycobacteria are acid fast
    - Does not provide species identification
  - Cannot be done without a culture

---

---

---

---

---

---

---

---

- Microscopy Results Guide Decisions**
- Clinical Management
    - Patient therapy may be initiated for TB
    - Changes in smear status important to monitor response to therapy
  - Public health interventions
    - Smear status and grade useful for identifying the most infectious cases
    - Contact investigation priority based on smear result
    - Decisions regarding respiratory isolation

---

---

---

---

---

---

---

---

### Direct Detection using PCR

- AKA: Nucleic Acid Amplification Testing (NAAT)

The diagram illustrates the PCR process. It starts with a sputum sample, followed by DNA extraction, and then amplification. The amplification step results in two possible outcomes: 'MTBC or MAC detected' and 'MTBC or MAC not detected'.

---

---

---

---

---

---

---

---

### PCR testing

- Detect *M. tuberculosis* complex (MTBC) and *M. avium* complex (MAC) directly from sputum sediment
- PCR test developed at WSLH (not-FDA cleared)
  - 1 other laboratory in the state is performing NAAT (Cepheid GeneXpert)
- Testing takes about 2 hours
- Unable to distinguish live and dead bacilli

---

---

---

---

---

---

---

---

### PCR testing (cont.)

- Automatically performed on new smear positive patients
- Fee-exempt testing for smear positive specimens and patients suspected of having active TB (approved by WI TB Program)
- Sensitivity
  - >95% for AFB smear-positive, culture-confirmed TB patients
  - 55-75% of AFB smear-negative, culture-confirmed TB patients

---

---

---

---

---

---

---

---